



User Guide

How to use the Breath Test system - in easy steps.

For nurses, lab technicians, doctors, and administrators.

This guide shows you each task, one step at a time. You do not need any computer knowledge to follow it. If you cannot see a button that this guide mentions, it simply means that job is not part of your role.

1. What is this system?

It is a program for hospitals that use QuinTron breath test machines. It helps you save a patient's breath test, turn the readings into a chart, and make a neat report that a doctor can check and sign.

What you can do with it

- Add a patient.
- Start a breath test (for example: SIBO Glucose, Lactose, or Fructose).
- Type in the breath readings taken during the test.
- See the results drawn as a chart, made for you.
- Let the doctor write the result and sign it.
- Save the report as a PDF, share it, or save the data as a spreadsheet.

Important: The system never decides the result on its own. It may show a few helpful hints, but the final result is always written and signed by a doctor.

2. Logging in

1. Open the web page your hospital gave you.
2. Type your email and password.
3. Click "Sign in".

Tip: If it will not let you in, check your email and password are spelled correctly. If it still will not work, ask your administrator.

3. Who can do what

Role	What they can do
Nurse / Lab Technician	Add patients, start tests, type in readings, and share reports.
Doctor	Check a test, write the result, and sign it. Can reopen a signed report if needed.
Administrator	Manage staff accounts, test types, hospitals, and view the activity list.
Specter Support	Help with setup. Can look at reports but cannot sign them.

4. Add a patient

1. On the home page, click "Patients".
2. Click "Register patient".
3. Fill in the boxes: patient ID (their hospital number), name, date of birth, gender, weight, hospital, and referring doctor.
4. Click "Register patient" to save.

Tip: Patient details are kept private and safe. Only the right staff can see them.

5. Start a breath test

1. Open the patient, then click "New breath test".
2. Pick the test type. The substrate and dose fill in for you, but you can change them if you need to.
3. Add the department, the dates, and the nurse or technician doing the test.
4. Write down any symptoms or notes from before the test.
5. Click "Create test". The test now shows as a "Draft".

What the test labels mean

- **Draft:** just started, no readings yet.
- **In Progress:** readings added, waiting for the doctor.
- **Finalized:** signed by the doctor and locked.

6. Type in the readings

1. Open the test and click "Enter samples".
2. Each row is one reading. Type the row number, the time in minutes, and the two gas numbers (H2 and CH4) from the machine.
3. The system adds H2 and CH4 together for you. You do not need to do any maths.
4. Fill in the other machine readings and any symptoms if you have them.
5. Click "+ Add sample" to add another row. Click "Save samples" when you are done.

If a reading was skipped

Tick the "Skip" box on that row and write the reason (for example, "patient stepped out"). A short reason is needed so nothing is left unexplained.

Tip: You can go back and change the readings at any time, right up until the doctor signs the report.

7. See the chart

As soon as you save some readings, the test page draws a chart for you. It shows how the two gases and their total change over time. The same chart appears in the report.

8. Writing the result and signing (doctor only)

1. Open the test and go to the doctor's review section.
2. The doctor writes the result and any advice for the patient.
3. Click "Save" to keep it for later, or "Sign and Finalize" to finish and lock the report.

Important: Once a report is signed, it cannot be changed. If a change is really needed, a doctor or administrator can "Reopen" it and give a reason. That reason is kept on record.

The system may show a few helpful hints, but they are only there to help. The doctor always makes the final decision.

9. Save, share, and print

Save the report as a PDF

1. Open the test.
2. Click "Download PDF". The full report opens, ready to print or save.

Share the report

- Click "Share" to copy a link or start an email.
- The person you send it to needs their own login to open it.

Save the data as a spreadsheet

- On a test, click "Export" to save that test's readings.
- On the tests list, click "Export" to save a summary of many tests.

Tip: These files open straight away in Microsoft Excel.

10. Administrator tasks

Only administrators can see these (and some of them, Specter Support).

Staff accounts

- Add a staff member with their name, email, role, and a first password.
- Change someone's role, or turn an account on or off.

Test types

- Add your own test types.
- Set the substrate and dose that fill in by default.

Hospitals and departments

- Add hospitals, and add departments inside each hospital.

Activity list

- See who did what and when: logins, changes, and sign-offs.

11. Safety and privacy

- Patient details are kept private and safe.
- Every important action is written down automatically.
- The system never decides a result by itself. The doctor is always in charge.

12. Need help?

If something does not work, or a button you expect is not there, ask your hospital administrator or Specter Support. Most of the time, a missing button just means that job is not part of your role.